

# Raw Data Report

May 2020



## Project Information

Client Name	MacroGen Oceania PL
Company / Institution	MacroGen Oceania PL
Order Number	HN00126866
Type of Read	Paired-end
Read Length	151
Number of Samples	18
Library Kit	TruSeq RNA Sample Prep Kit v2
Library Protocol	TruSeq RNA Sample Preparation v2 Guide, Part # 15026495 Rev. F
Type of Sequencer	Illumina platform

# Table of Contents

---

<b>Project Information</b>	<b>02</b>
<b>1. Data Download Information</b>	<b>04</b>
1. 1. Raw Data and Analysis Results	04
<b>2. Experimental Methods and Workflow</b>	<b>06</b>
2. 1. Experiment Overview	06
2. 2. Generation of Raw Data	07
<b>3. Summary of Produced Data</b>	<b>08</b>
3. 1. Raw Data Statistics	08
3. 2. Total Read Bases	09
3. 3. Total Reads	10
3. 4. GC/AT Content	11
3. 5. Q20/Q30 (%)	12
<b>4. Appendix</b>	<b>13</b>
4. 1. FAQ	13
4. 2. FASTQ File	13
4. 3. Phred Quality Score Chart	13

# 1. Data Download Information

## 1. 1. Raw Data and Analysis Results

Download link	File size	md5sum
<a href="#">P1A2_1.fastq.gz</a>	1.9G	e1b119ef028b371b187a0cd0ae1685f2
<a href="#">P1A2_2.fastq.gz</a>	1.9G	714ca007d6d3012811ba6c6f99848a81
<a href="#">P1A3_1.fastq.gz</a>	1.9G	0c839bca5e505bb2eb6c909c0d0e2c9a
<a href="#">P1A3_2.fastq.gz</a>	1.9G	093c76f39b430a7def08141f8a69521e
<a href="#">P1B2_1.fastq.gz</a>	1.9G	328a97fbe7143beade1f4d68bcc4d86
<a href="#">P1B2_2.fastq.gz</a>	2.0G	22e02a9121f885666e1f4cfe74d5112
<a href="#">P1B3_1.fastq.gz</a>	1.6G	5efbb4b7ec12d811724175d8fe486496
<a href="#">P1B3_2.fastq.gz</a>	1.6G	4481c63d454a33c8725a3aa3cace70da
<a href="#">P1C2_1.fastq.gz</a>	1.9G	05149fb45a90f133708f3f4c8c918d69
<a href="#">P1C2_2.fastq.gz</a>	1.9G	59b540ffdd73279d06748f830fe021fd
<a href="#">P1C3_1.fastq.gz</a>	1.6G	d0ce85227a8dd3cc076e65395516789a
<a href="#">P1C3_2.fastq.gz</a>	1.7G	2ba4e5a12d7dc9c4e94692578d83246b
<a href="#">P2A2_1.fastq.gz</a>	1.9G	ba334190b2f37bda5364b23da6019d66
<a href="#">P2A2_2.fastq.gz</a>	1.9G	c1d5bad7e7c01f0c0bc3d56026feb90
<a href="#">P2A3_1.fastq.gz</a>	1.6G	05d912c4ab91e8f925b76f60585ea9cb
<a href="#">P2A3_2.fastq.gz</a>	1.7G	9a87fc06635e396515ad7f214c643834
<a href="#">P2B2_1.fastq.gz</a>	1.9G	6ab1fa1c19779033e2c076cffcbf2541
<a href="#">P2B2_2.fastq.gz</a>	1.9G	20b9d9cb92cb3ab87492212469db755b
<a href="#">P2B3_1.fastq.gz</a>	1.9G	6c73b72dc427dbefa22cea8d8332c384
<a href="#">P2B3_2.fastq.gz</a>	2.0G	3ee03d49ed4daa294942d0e3fa68d8e1
<a href="#">P2C2_1.fastq.gz</a>	1.9G	51f6fabb9ad6bcd5c55eb356aca1b421
<a href="#">P2C2_2.fastq.gz</a>	1.9G	2d67871cb7626c1132fe500f72245468
<a href="#">P2C3_1.fastq.gz</a>	1.8G	7ed061f4a1001fc330302c6efbc8b570
<a href="#">P2C3_2.fastq.gz</a>	1.9G	3a3ef8afd308331e9861874992ee524c
<a href="#">P3A2_1.fastq.gz</a>	1.9G	ea3a0ec2be04e5ce4824bae525c46d32
<a href="#">P3A2_2.fastq.gz</a>	1.9G	891a3e9561ee97bc1d5a3c46d0a5e01e
<a href="#">P3A3_1.fastq.gz</a>	2.1G	40eba46ba3cbe4dd02a9ccd86d282b2b
<a href="#">P3A3_2.fastq.gz</a>	2.2G	1ca639b80c39d9b97a23f4b359cff6fe
<a href="#">P3B2_1.fastq.gz</a>	1.8G	535672cdaed667246e6f469d7fbee66a
<a href="#">P3B2_2.fastq.gz</a>	1.9G	a3340bf332b5231c7c838aca7bfa9bb7
<a href="#">P3B3_1.fastq.gz</a>	1.9G	406b9a83d65ceaa3bd1e0f19452b4683
<a href="#">P3B3_2.fastq.gz</a>	2.0G	8b677a384e847a4bf9f549c4610ee32e
<a href="#">P3C2_1.fastq.gz</a>	1.5G	d586532625d0a0a4d1d13792c62af609

<a href="#">P3C2_2.fastq.gz</a>	1.6G	8e20bfbcc96d3d122697e8d0708ad9fa
<a href="#">P3C3_1.fastq.gz</a>	1.9G	9d1b89ac658b41b1d728b6a3e4443894
<a href="#">P3C3_2.fastq.gz</a>	1.9G	f392db46588deb8ee8481735fd296bc5

- **fastq.gz** : This is a zip file of raw data used in analysis.
- **md5sum** : In order to verify the integrity of files, md5sum is used. If the values of md5sum are the same, there is no forgery, modification or omission.

**Your data will be retained in our server for 3 months. Should you wish to extend the retention period, please contact us.**

## 2. Experimental Methods and Workflow

### 2. 1. Experiment Overview



Fig1. Experiment overview

The Illumina NGS workflow includes 4 basic steps :

#### 1) Sample Preparation

For library construction, DNA/RNA is extracted from a sample. After performing quality control (QC), qualified samples proceed to library construction.

#### 2) Library Construction

The sequencing library is prepared by random fragmentation of the DNA or cDNA sample, followed by 5' and 3' adapter ligation. Alternatively, "tagmentation" combines the fragmentation and ligation reactions into a single step that greatly increases the efficiency of the library preparation process. Adapter-ligated fragments are then PCR amplified and gel purified.

#### 3) Sequencing

For cluster generation, the library is loaded into a flow cell where fragments are captured on a lawn of surface-bound oligos complementary to the library adapters. Each fragment is then amplified into distinct, clonal clusters through bridge amplification. When cluster generation is complete, the templates are ready for sequencing.

Illumina SBS technology utilizes a proprietary reversible terminator-based method that detects single bases as they are incorporated into DNA template strands. As all 4 reversible, terminator-bound dNTPs are present during each sequencing cycle, natural competition minimizes incorporation bias and greatly reduces raw error rates compared to other technologies. The result is highly accurate base-by-base sequencing that virtually eliminates sequence-context-specific errors, even within repetitive sequence regions and homopolymers.

#### 4) Raw data

Sequencing data is converted into raw data for the analysis.

## 2. 2. Generation of Raw Data

The Illumina sequencer generates raw images utilizing sequencing control software for system control and base calling through an integrated primary analysis software called RTA (Real Time Analysis). The BCL (base calls) binary is converted into FASTQ utilizing illumina package bcl2fastq. Adapters are not trimmed away from the reads.

## 3. Summary of Produced Data

### 3. 1. Raw Data Statistics

The total number of bases, reads, GC (%), Q20 (%), and Q30 (%) are calculated for the 18 samples. For example, in P1A2, 59,887,438 reads are produced, and total read bases are 9.0G bp. The GC content (%) is 50.35% and Q30 is 94.8%.

Table 1. Raw data Stats (maximum 20 samples)

Sample ID	Total read bases (bp)	Total reads	GC(%)	AT(%)	Q20(%)	Q30(%)
P1A2	9,043,003,138	59,887,438	50.35	49.65	98.28	94.8
P1A3	8,976,955,436	59,450,036	50.4	49.6	98.35	94.95
P1B2	8,134,225,946	53,869,046	50.02	49.98	98.19	94.55
P1B3	6,844,404,784	45,327,184	50.29	49.71	98.3	94.85
P1C2	9,008,294,580	59,657,580	50.44	49.56	98.18	94.59
P1C3	7,098,403,696	47,009,296	50.8	49.2	98.31	94.91
P2A2	9,018,764,014	59,726,914	49.75	50.25	98.26	94.8
P2A3	7,072,656,686	46,838,786	50.27	49.73	98.3	94.84
P2B2	8,965,561,882	59,374,582	50.55	49.45	98.23	94.7
P2B3	8,242,593,512	54,586,712	49.85	50.15	98.25	94.72
P2C2	8,952,539,944	59,288,344	50.13	49.87	98.26	94.75
P2C3	7,860,717,230	52,057,730	50.37	49.63	98.28	94.79
P3A2	8,972,362,016	59,419,616	49.78	50.22	98.23	94.69
P3A3	9,032,326,532	59,816,732	49.84	50.16	98.25	94.76
P3B2	7,796,249,290	51,630,790	50.33	49.67	98.31	94.88
P3B3	8,174,337,888	54,134,688	50.58	49.42	98.27	94.75
P3C2	6,540,168,172	43,312,372	50.79	49.21	98.32	94.9
P3C3	8,978,070,118	59,457,418	49.47	50.53	98.3	94.85

- Sample ID : Sample name.
- Total read bases : Total number of bases sequenced.
- Total reads : Total number of reads. For Illumina paired-end sequencing, this value refers to the sum of read 1 and read 2.
- GC(%) : GC content.
- AT(%) : AT content.
- Q20(%) : Ratio of bases that have phred quality score of over 20.
- Q30(%) : Ratio of bases that have phred quality score of over 30.



## 3. 2. Total Read Bases




Figure 2. Throughput of Raw data

### 3. 3. Total Reads




Figure 3. Total read count of Raw data

### 3. 4. GC/AT Content




Figure 4. GC/AT Content of Raw data

### 3. 5. Q20/Q30 (%)




Figure 5. Q20/Q30 scores of Raw data

## 4. Appendix

### 4. 1. FAQ

**Q:** I want to see the produced data. How can I open the files?

**A:** As the large size zip files provided by our company are hard to process in the Windows environment, we highly recommend using Linux environment for a smoother operation.

### 4. 2. FASTQ File

Example of FASTQ

```
@HISEQ-MFG:501:HB0TFADXX:1:1101:1247:2183 1:N:0:
CTCAGCTAAATACTTTGACACCNGTANNANNNNNNNNNNTNNNNNNNNNNNN
+
@@@BDDDDHHHHFHIIIIIII#3AC#####
```

FASTQ file is composed of four lines.

Line 1 : ID line includes information such as flow cell lane information.

Line 2 : Sequences line.

Line 3 : Separator line (+ mark).

Line 4 : Quality values line about sequences.

### 4. 3. Phred Quality Score Chart

Phred quality score numerically expresses the accuracy of each nucleotide. Higher Q number signifies higher accuracy. For example, if Phred assigns a quality score of 30 to a base, the chances of having base call error are 1 in 1000.

Phred Quality Score Q is calculated with  $-10\log_{10}P$ , where P is probability of erroneous base call.

Quality of phred score	Probability of incorrect base call	Base call accuracy	Characters
10	1 in 10	90%	!"#\$%&'()*+,-./012345
20	1 in 100	99%	6789;:h=i?
30	1 in 1000	99.9%	@ABCDEFGHIJ
40	1 in 10000	99.99%	

- Encoding : Sanger Quality (ASCII Character Code=Phred Quality Value + 33)



## HEADQUARTER

### Macrogen, Inc.

#### Laboratory, IT and Business Headquarter & Support Center

[08511] 1001, 10F, 254, Beotkkot-ro,  
Geumcheon-gu, Seoul, Republic of Korea  
(Gasan-dong, World Meridian 1)  
Tel: +82-2-2180-7000  
Email1: ngs@macrogen.com(Overseas)  
Email2: ngskr@macrogen.com  
(Republic of Korea)  
Web: www.macrogen.com  
LIMS: dna.macrogen.com

## SUBSIDIARY

### Macrogen Europe

#### Laboratory, Business & Support Center

Meibergdreef 31, 1105 AZ, Amsterdam,  
the Netherlands  
Tel: +31-20-333-7563  
Email: ngs@macrogen.eu

### Macrogen Singapore

#### Laboratory, Business & Support Center

3 Biopolis Drive #05-18, Synapse,  
Singapore 138623  
Tel: +65-6339-0927  
Email: info-sg@macrogen.com

## BRANCH

### Macrogen Spain

#### Laboratory, Business & Support Center

Av. Sur del Aeropuerto de Barajas,  
28. Office B-2, 28042 Madrid, Spain  
Tel: +34-911-138-378  
Email: info-spain@macrogen.com

### Psomagen (Macrogen USA)

#### Laboratory, Business & Support Center

1330 Piccard Drive, Suite 103, Rockville,  
MD 20850, United States  
Tel: +1-301-251-1007  
Email: inquiry@psomagen.com

### Macrogen Japan

#### Laboratory, Business & Support Center

3F Kyoto University International Science  
Innovation Bldg.  
36-1 Yoshida-honmachi, Sakyo-ku,  
Kyoto 606-8501 JAPAN  
Tel: +81-75-746-2773  
Email: customer@macrogen-japan.co.jp